

2021年12月

スフィンゴイド塩基

長良サイエンス株式会社

〒501-1121 岐阜市古市場 840


Tel 058-234-4257, Fax 058-234-4724

E-mail: nagara@nsgifu.jp, <http://www.nsgifu.jp>


Code No.	Storage Conditions	Product Name · Purity · Molecular Formula = Molecular Weight · [CAS Registry Number]	Quantity · Price (JPY)
NS440101	-20 °C	<i>trans</i> -4- <i>cis</i> -8-Sphingadienine ≥98% (HPLC)	1mg 40,000
NS440102		C ₁₈ H ₃₅ NO ₂ = 297.48 [41679-33-8]	5mg 120,000
NS440801	-20 °C	<i>trans</i> -4- <i>trans</i> -8-Sphingadienine ≥98% (HPLC)	1mg 50,000 ←新製品
		C ₁₈ H ₃₅ NO ₂ = 297.48 [40878-79-3]	
NS440201	-20 °C	4-Hydroxy- <i>cis</i> -8-sphingenine ≥98% (HPLC)	1mg 65,000
NS440202		C ₁₈ H ₃₇ NO ₃ = 315.49 [218786-92-6]	5mg 200,000
NS440901	-20 °C	9-Methyl- <i>trans</i> -4- <i>trans</i> -8-sphingadienine ≥98% (HPLC)	1mg 60,000 ←新製品
		C ₁₉ H ₃₇ NO ₂ = 311.50 [160226-28-8]	
NS440302	-20 °C	C ₁₂ -Sphingosine ≥99% (HPLC)	10mg 36,000
		C ₁₂ H ₂₅ NO ₂ = 215.33 [128427-86-1]	
NS440402	-20 °C	C ₁₄ -Sphingosine ≥99% (HPLC)	10mg 36,000
		C ₁₄ H ₂₉ NO ₂ = 243.39 [24558-60-9]	
NS440502	-20 °C	C ₁₆ -Sphingosine ≥99% (HPLC)	10mg 36,000
		C ₁₆ H ₃₃ NO ₂ = 271.44 [6982-09-8]	
NS440602	-20 °C	C ₁₈ -Sphingosine ≥99% (HPLC)	25mg 14,000
		C ₁₈ H ₃₇ NO ₂ = 299.49 [123-78-4]	
NS440702	-20 °C	C ₂₀ -Sphingosine ≥99% (HPLC)	10mg 36,000
		C ₂₀ H ₄₁ NO ₂ = 327.55 [6918-49-6]	




trans-4-*cis*-8-Sphingadienine




trans-4-*trans*-8-Sphingadienine



4-Hydroxy-*cis*-8-sphingenine



9-Methyl-*trans*-4-*trans*-8-sphingadienine



R=C₇H₁₅, C₉H₁₉, C₁₁H₂₃, C₁₃H₂₇, C₁₅H₃₁